


Dual-Field Microfluidic Chip for Label-Free CTC Isolation

Integrating Inertial Focusing and Acoustophoresis

Carl Chen

Culver Academies

Challenges in CTC Isolation

A diagram illustrating the challenges in circulating tumor cell (CTC) isolation. On the left, a large, irregular cluster of brown, polygonal cells represents a tumor. A line points from the text 'Tumor Cells' to this cluster. On the right, a horizontal red band represents a blood vessel. Inside the vessel, several blue, elongated cells with dark blue nuclei are shown, representing CTCs. A line points from the text 'CTC' to one of these cells. The background is a dark grey gradient.

- CTCs are extremely rare (~1–10 cells/mL)
- Current methods rely on antibody labeling (costly, marker-dependent)
- Need for label-free, high-throughput separation

Inertial separation





Image made by Student

Acoustic separation



- RBC (red blood cell): 6-8 μm
- WBC (white blood cell): 10-12 μm
- CTC (circulating tumor cell): 10-20 μm

- Combine inertial microfluidics + acoustic separation
- Sequential, dual-field design
- Add anti-adhesion strategy
- Sequential physical separation improves both efficiency and stability.

Engineering Approach

Chip Design Overview

- Spiral channel (inertial)
- Straight channel (acoustic)
- PZT transducer
- Outlet




Image made by Student

$$Re = \frac{\rho U D_h}{\mu}$$


Stage 1: Inertial Focusing

- Dean flow in spiral channel
- Size-dependent lateral migration
- Larger particles move toward equilibrium positions
- Lift force \propto particle size⁴




$$F_L \propto \rho U^2 a^4 / D_h^2$$

Stage 2: Acoustic Separation




- Standing wave field
- Acoustic radiation force
- Larger particles migrate to pressure nodes



Anti-Fouling Strategy

- SiO₂ micropillar structures
- PEG coating
- Reduces particle-wall adhesion
- Improves long-term stability



Experimental Results

- Inertial:

- 96% enrichment
- Stable flow behavior

- Acoustic:

- OFF → mixed distribution
- ON → clear focusing




Image made by Student




Image made by Student





Image made by Student

Limitations & Future Work

- No live-cell validation yet
- Need statistical analysis (larger sample size)
- Optimize acoustic parameters
- Future: clinical testing




Image made by Student

Conclusion & Impact

- Label-free separation platform
- Scalable, cost-effective
- Potential for precision medicine




Image made by Student